



Minimum System Requirements

Use the guidelines in this document as a guide to needed equipment before making any purchases.

Operating Systems

Windows 10 (64 Bit)

Additional System Requirements

Windows.NET Framework 3.5 SP1 and .NET 4.0

Display

1280 x 1024 resolution or better

24-bit color/32-bit color depth recommended

Open GL 1.1 compliance

CPU

Intel or AMD compatible 2.6GHz or better/Dual or Quad core 2.6 GHz or better CPUs recommended

Memory

64-bit systems: 16GB recommended minimum

Free Disk Space

10GB or more

GPU

NVIDIA GPU less than 2 years old with CUDA. The NVIDIA driver must be > 418.96